



# **Dil Membrane Staining Kit**

## **User Manual**

**Catalog # CRG1067**

Highly orange fluorescent membrane probe for labeling live and fixed cells


**For research use only. Not for diagnostic or therapeutic procedures.**

I. INTRODUCTION.....	2
II. MATERIALS REQUIRED BUT NOT PROVIDED.....	3
III. PROCEDURAL GUIDELINES.....	3
IV. WORKING SOLUTION PREPARATION.....	3
V. ASSAY PROCEDURE.....	4
VI. TECHNICAL SUPPORT.....	5
VII. NOTES.....	5

## I. INTRODUCTION

Dil (1,1'-dioctadecyl-3,3'-tetramethylindocarbocyanine perchlorate) has the molecular formula of  $C_{59}H_{97}ClN_2O_4$  with a molecular weight of 933.88. Its CAS Number is 41085-99-8. Dil is one of the most commonly used cell membrane fluorescent probes, exhibiting orange-red fluorescence when bound to cell membrane. Dil is a lipophilic carbocyanine dye with a long lipophilic hydrocarbon chain, enabling its lateral diffusion in cell membrane.

Only Dil bound to cell membrane can be excited to emit strong orange-red fluorescence. Please refer to Figure 1 for the excitation and emission spectra of Dil after binding to phosphate bilayer of cell membrane.



Dil is widely used as a tracer or long-term tracer for positive or retrograde, live or fixed cells or tissues. It has low toxicity and usually does not affect cell viability, thereby Dil-labeled neuron cells can survive up to 4 weeks in vitro and up to a year in vivo. The migration rate of Dil is 0.2-0.6 mm/day in neuron cell membranes that have undergone fixation and 6 mm/day in living neuron cell membranes.

The Staining Enhancer provided in this kit enables faster staining of cell membranes, brighter fluorescence, and lower staining background.

Beside the fluorescent labeling of cell membranes, this product can also be used to detect cell fusion and adhesion, cell migration during development or transplantation, lipid diffusion across cell membranes by FRAP (Fluorescence Recovery After Photobleaching), cytotoxicity and lipoproteins, etc.

## II. KIT COMPONENTS

Component	50 Assays	300 Assays	Storage
Dil Probe (100X)	50 $\mu$ l x 1	300 $\mu$ l x 1	4 °C
Dilution Buffer	5 ml x 1	30 ml x 1	4 °C
Technical Manual	1 Manual	1 Manual	

## III. STORAGE AND STABILITY

Shipped at 4°C. Store at -20°C and protect from light for 12 months.

## IV. WORKING SOLUTION PREPARATION

Working Solution: Dilute 1  $\mu$ l Dil Probe (100X) in 100  $\mu$ l Dilution Buffer.

*Note:* The dilution ratio can be adjusted appropriately according to the experimental effect.

Ex/Em = 549/565 nm


## V. ASSAY PROCEDURE

### 1. For adherent cells staining

- (1) Adherent cells were cultured on a sterile cover slide.
- (2) Remove the cover glass from the medium, absorb excess liquid but keep the surface moist.
- (3) Add 100  $\mu$ L of Working Solution to one corner of the cover glass and gently shake to evenly cover all cells.
- (4) Cells were incubated at 37 °C for 2-20 mins. The reaction time can be optimized to obtain uniform labeling effect.
- (5) Discard Working Solution, wash the glass with PBS for 2 to 3 times.

### 2. For suspension cells staining

- (1) Adding an appropriate volume of Working Solution to re-suspension cells, the density of the cells is  $1 \times 10^6$  /mL.
- (2) The cells were incubated at 37°C for 2-20 min. The reaction time can be optimized to obtain uniform labeling effect.
- (3) After incubation, centrifuge at 1000-1500 rpm for 5 mins. Discard the supernatant and slowly add the growth medium again to resuspend the cells.
- (4) Repeat step (3) more than twice.



Immunofluorescent analysis of Dil staining in PC3 cells. Formalin-fixed cells were permeabilized with 0.1% Triton X-100 in TBS for 5-10 minutes and blocked with 3% BSA-PBS for 30 minutes at room temperature. Cells were probed with the primary antibody in 3% BSA-PBS and incubated overnight at 4 °C in a humidified chamber. Cells were washed with PBST and incubated with Dil (orange) at room temperature in the dark. Phalloidin - AF488 was used to stain Actin filaments (green). DAPI was used to stain the cell nuclei (blue).

## **VI. TECHNICAL SUPPORT**

For troubleshooting, information or assistance, please go online to [www.cohesionbio.com](http://www.cohesionbio.com) or contact us at [techsupport@cohesionbio.com](mailto:techsupport@cohesionbio.com)

## **VII. NOTES**