

# Product Data Sheet

## CHST2 Blocking Peptide

| Catalog #                 | Source                                                                                                                                                                                          | Reactivity | Applications |
|---------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|--------------|
| CBP2391                   | Synthetic                                                                                                                                                                                       | H          | BL           |
| <b>Description</b>        | The peptide is used to block Anti-CHST2 Antibody (#CPA2391) reactivity.                                                                                                                         |            |              |
| <b>Form</b>               | Lyophilized powder                                                                                                                                                                              |            |              |
| <b>Gene Symbol</b>        | CHST2                                                                                                                                                                                           |            |              |
| <b>Alternative Names</b>  | GN6ST; Carbohydrate sulfotransferase 2;<br>Galactose/N-acetylglucosamine/N-acetylglucosamine 6-O-sulfotransferase 2; GST-2;<br>N-acetylglucosamine 6-O-sulfotransferase 1; GlcNAc6ST-1; Gn6ST-1 |            |              |
| <b>Entrez Gene</b>        | 9435 (Human)                                                                                                                                                                                    |            |              |
| <b>SwissProt</b>          | Q9NYF8 (Human)                                                                                                                                                                                  |            |              |
| <b>Purity</b>             | >85%                                                                                                                                                                                            |            |              |
| <b>Quality Control</b>    | The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.                                                                                                       |            |              |
| <b>Directions for Use</b> | Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.               |            |              |
| <b>Storage/Stability</b>  | Shipped at 4°C. Store at -20°C for one year.                                                                                                                                                    |            |              |

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

### COHESION BIOSCIENCES LIMITED

**WEB**  
www.cohesionbio.com

**ORDER**  
order@cohesionbio.com

**SUPPORT**  
techsupport@cohesionbio.com

**CUSTOM**  
custom@cohesionbio.com