


## Angiotensin I [Des-Asp1], Human

Catalog #	Source	Reactivity	Applications
CCP1072	Synthetic		
<b>Description</b>	Peptide to Angiotensin I [Des-Asp1], Human		
<b>Biological Description</b>			
<b>Form</b>	Lyophilized powder		
<b>CAS Number</b>	56317-01-2		
<b>Molecular Formula</b>	C <sub>58</sub> H <sub>84</sub> N <sub>16</sub> O <sub>11</sub>		
<b>Molecular Weight</b>	1181.42		
<b>Purity</b>	> 95%		
<b>Chemical Structure</b>	H - Arg - Val - Tyr - Ile - His - Pro - Phe - His - Leu - OH		
<b>Storage/Stability</b>	Shipped at 4°C. Store at -20°C for one year.		

### Structure



**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

### COHESION BIOSCIENCES LIMITED

**WEB**  
www.cohesionbio.com

**ORDER**  
order@cohesionbio.com

**SUPPORT**  
techsupport@cohesionbio.com

**CUSTOM**  
custom@cohesionbio.com