


Atriopeptin II, Mouse, Rat, Rabbit

Catalog #	Source	Reactivity	Applications
CCP1091	Synthetic		
Description	Peptide to Atriopeptin II, Mouse, Rat, Rabbit		
Biological Description			
Form	Lyophilized powder		
CAS Number	89139-54-8		
Molecular Formula	C98H156N34O32S2		
Molecular Weight	2386.67		
Purity	> 95%		
Chemical Structure	H - Ser - Ser - Cys - Phe - Gly - Gly - Arg - Ile - Asp - Arg - Ile - Gly - Ala - Gln - Ser - Gly - Leu - Gly - Cys - Asn - Ser - Phe - Arg - OH (Disulfide bridge Cys7 - Cys23)		
Storage/Stability	Shipped at 4°C. Store at -20°C for one year.		

Structure



Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com