

Helospectin II

Catalog #	Source	Reactivity	Applications
CCP1255	Synthetic		
Description	Peptide to Helospectin II		
Biological Description			
Form	Lyophilized powder		
CAS Number			
Molecular Formula	C180H288N46O57		
Molecular Weight	4008.58		
Purity	> 95%		
Chemical Structure	H - His - Ser - Asp - Ala - Thr - Phe - Thr - Ala - Glu - Tyr - Ser - Lys - Leu - Leu - Ala - Lys - Leu - Ala - Leu - Gln - Lys - Tyr - Leu - Glu - Ser - Ile - Leu - Gly - Ser - Ser - Thr - Ser - Pro - Arg - Pro - Pro - Ser - OH		
Storage/Stability	Shipped at 4°C. Store at -20°C for one year.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com