

## Melanostatin, Frog

Catalog #	Source	Reactivity	Applications
CCP1314	Synthetic		
<b>Description</b>	Peptide to Melanostatin, Frog		
<b>Biological Description</b>			
<b>Form</b>	Lyophilized powder		
<b>CAS Number</b>			
<b>Molecular Formula</b>	C189H285N53O57S		
<b>Molecular Weight</b>	4243.76		
<b>Purity</b>	> 95%		
<b>Chemical Structure</b>	H - Tyr - Pro - Ser - Lys - Pro - Asp - Asn - Pro - Gly - Glu - Asp - Ala - Pro - Ala - Glu - Asp - Met - Ala - Lys - Tyr - Tyr - Ser - Ala - Leu - Arg - His - Tyr - Ile - Asn - Leu - Ile - Thr - Arg - Gln - Arg - Tyr - NH2		
<b>Storage/Stability</b>	Shipped at 4°C. Store at -20°C for one year.		

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

### COHESION BIOSCIENCES LIMITED

**WEB**  
[www.cohesionbio.com](http://www.cohesionbio.com)

**ORDER**  
[order@cohesionbio.com](mailto:order@cohesionbio.com)

**SUPPORT**  
[techsupport@cohesionbio.com](mailto:techsupport@cohesionbio.com)

**CUSTOM**  
[custom@cohesionbio.com](mailto:custom@cohesionbio.com)