

**Orexin B, Canine**

Catalog #	Source	Reactivity	Applications
CCP1359	Synthetic		
<b>Description</b>	Peptide to Orexin B, Canine		
<b>Biological Description</b>			
<b>Form</b>	Lyophilized powder		
<b>CAS Number</b>			
<b>Molecular Formula</b>	C125H214N44O34S		
<b>Molecular Weight</b>	2909.44		
<b>Purity</b>	> 95%		
<b>Chemical Structure</b>	H - Arg - Pro - Gly - Pro - Pro - Gly - Leu - Gln - Gly - Arg - Leu - Gln - Arg - Leu - Leu - Gln - Ala - Ser - Gly - Asn - His - Ala - Ala - Gly - Ile - Leu - Thr - Met - NH <sub>2</sub>		
<b>Storage/Stability</b>	Shipped at 4°C. Store at -20°C for one year.		

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

**COHESION BIOSCIENCES LIMITED****WEB**[www.cohesionbio.com](http://www.cohesionbio.com)**ORDER**[order@cohesionbio.com](mailto:order@cohesionbio.com)**SUPPORT**[techsupport@cohesionbio.com](mailto:techsupport@cohesionbio.com)**CUSTOM**[custom@cohesionbio.com](mailto:custom@cohesionbio.com)