

Orexin B, Mouse, Rat

Catalog #	Source	Reactivity	Applications
CCP1361	Synthetic		
Description	Peptide to Orexin B, Mouse, Rat		
Biological Description			
Form	Lyophilized powder		
CAS Number	202801-92-1		
Molecular Formula	C126H215N45O34S		
Molecular Weight	2936.46		
Purity	> 95%		
Chemical Structure	H - Arg - Pro - Gly - Pro - Pro - Gly - Leu - Gln - Gly - Arg - Leu - Gln - Arg - Leu - Leu - Gln - Ala - Asn - Gly - Asn - His - Ala - Ala - Gly - Ile - Leu - Thr - Met - NH ₂		
Storage/Stability	Shipped at 4°C. Store at -20°C for one year.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com