


Secretin, Human

Catalog #	Source	Reactivity	Applications
CCP1416	Synthetic		
Description	Peptide to Secretin, Human		
Biological Description	Gastrointestinal hormone which stimulates pancreatic secretion.		
Form	Lyophilized powder		
CAS Number	108153-74-8		
Molecular Formula	C130H220N44O40		
Molecular Weight	3039.47		
Purity	> 95%		
Chemical Structure	H - His - Ser - Asp - Gly - Thr - Phe - Thr - Ser - Glu - Leu - Ser - Arg - Leu - Arg - Glu - Gly - Ala - Arg - Leu - Gln - Arg - Leu - Leu - Gln - Gly - Leu - Val - NH ₂		
Storage/Stability	Shipped at 4°C. Store at -20°C for one year.		
Structure			



Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com