


[Nle13]-Motillin

Catalog #	Source	Reactivity	Applications
CCP1497	Synthetic		
Description	Peptide to [Nle13]-Motillin		
Biological Description			
Form	Lyophilized powder		
CAS Number			
Molecular Formula	C121H190N34O35		
Molecular Weight	2681.07		
Purity	> 95%		
Chemical Structure	Phe - Val - Pro - Ile - Phe - Thr - Tyr - Gly - Glu - Leu - Gln - Arg - Nle - Gln - Glu - Lys - Glu - Arg - Asn - Lys - Gly - Gln		
Storage/Stability	Shipped at 4°C. Store at -20°C for one year.		

Structure



Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com