

[Lys8]-Conopressin S

Catalog #	Source	Reactivity	Applications
CCP1605	Synthetic		

Description Peptide to [Lys8]-Conopressin S

Biological Description

Form Lyophilized powder

CAS Number 130836-26-9

Molecular Formula C₄₁H₇₃N₁₅O₁₀S₂


Molecular Weight 1000.26

Purity > 95%

Chemical Structure Cys - Ile - Ile - Arg - Asn - Cys - Pro - Lys - Gly - NH₂ (Disulfide bridge Cys1 - Cys6)

Storage/Stability Shipped at 4°C. Store at -20°C for one year.

Structure



Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com