


[Val5]-Angiotensin I (Human)

Catalog #	Source	Reactivity	Applications
CCP1671	Synthetic		
Description	Peptide to [Val5]-Angiotensin I (Human)		
Biological Description			
Form	Lyophilized powder		
CAS Number	484-43-5		
Molecular Formula	C ₆₁ H ₈₇ N ₁₇ O ₁₄		
Molecular Weight	1282.50		
Purity	> 95%		
Chemical Structure	Asp - Arg - Val - Tyr - Val - His - Pro - Phe - His - Leu		
Storage/Stability	Shipped at 4°C. Store at -20°C for one year.		

Structure



Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com