

Ac-Adhesin (1025-1044) amide

Catalog #	Source	Reactivity	Applications
CCP1740	Synthetic		
Description	Peptide to Ac-Adhesin (1025-1044) amide		
Biological Description			
Form	Lyophilized powder		
CAS Number	320350-56-9		
Molecular Formula	C97H160N26O32		
Molecular Weight	2202.51		
Purity	> 95%		
Chemical Structure	Ac - Gln - Leu - Lys - Thr - Ala - Asp - Leu - Pro - Ala - Gly - Arg - Asp - Glu - Thr - Thr - Ser - Phe - Val - Leu - Val - NH ₂		
Storage/Stability	Shipped at 4°C. Store at -20°C for one year.		
Structure			

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com