

## Mouse Anti-Monkey IgG (H&L) FITC

| Catalog #                | Source   | Reactivity | Applications |
|--------------------------|--|------------|--------------|
| CSA2339                  | Mouse  | Mk         | IF, FC       |
| <b>Description</b>       | Mouse Monoclonal Secondary Antibody to Monkey IgG (H&L) FITC   |            |              |
| <b>Immunogen</b>         | Monkey IgG   |            |              |
| <b>Purification</b>      | The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads.   |            |              |
| <b>Specificity</b>       | By immunoelectrophoresis and ELISA this antibody reacts specifically with Monkey IgG. Cross reactivity with human IgG. No antibody was detected against non immunoglobulin serum proteins. |            |              |
| <b>Clonality</b>         | Monoclonal   |            |              |
| <b>Conjugation</b>       | FITC   |            |              |
| <b>Form</b>              | Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.05% Sodium Azide   |            |              |
| <b>Dilution</b>          | IF (1/250 - 1/2000), FC (1/250 - 1/2000)   |            |              |
| <b>Storage/Stability</b> | Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.   |            |              |

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

### COHESION BIOSCIENCES LIMITED

**WEB**  
[www.cohesionbio.com](http://www.cohesionbio.com)

**ORDER**  
[order@cohesionbio.com](mailto:order@cohesionbio.com)

**SUPPORT**  
[techsupport@cohesionbio.com](mailto:techsupport@cohesionbio.com)

**CUSTOM**  
[custom@cohesionbio.com](mailto:custom@cohesionbio.com)