

Product Data Sheet

Anti-SAA Antibody (Capture)

Catalog #	Source	Reactivity	Applications
CAP1011-C	Mouse	H	ELISA, ICA, CLIA
Description	Mouse monoclonal capture antibody to SAA		
Immunogen	Recombinant protein of human SAA		
Purification	The antibody was purified by Protein A/G affinity chromatography.		
Specificity	Recognizes native and recombinant human SAA protein.		
Clonality	Monoclonal		
Conjugation			
Form	IgG1, liquid in PBS, pH 7.4.		
Dilution			
Gene Symbol	SAA1		
Alternative Names	Serum amyloid A-1 protein; SAA; Serum amyloid A-2 protein; SAA2		
Entrez Gene	6288		
SwissProt	P0DJ18		
Directions for Use	Recommended pair for sandwich immunoassay - Detection: CAP1011-D		
Storage/Stability	Shipped at 4°C. Upon delivery aliquot and store at -20°C for 2 years. Avoid freeze/thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com