

Product Data Sheet

Anti-Pepsinogen II Antibody (Detection)

Catalog #	Source	Reactivity	Applications
CAP1017-D	Mouse	H	ELISA, ICA, CLIA
Description	Mouse monoclonal detection antibody to Pepsinogen II		
Immunogen	Recombinant protein of human Pepsinogen II		
Purification	The antibody was purified by Protein A/G affinity chromatography.		
Specificity	Recognizes native and recombinant human Pepsinogen II protein.		
Clonality	Monoclonal		
Conjugation			
Form	IgG1, liquid in PBS, pH 7.4.		
Dilution			
Gene Symbol	PGC		
Alternative Names	Gastricsin; Pepsinogen C		
Entrez Gene	5225		
SwissProt	P20142		
Directions for Use	Recommended pair for sandwich immunoassay - Capture: CAP1017-C		
Storage/Stability	Shipped at 4°C. Upon delivery aliquot and store at -20°C for 2 years. Avoid freeze/thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com