

Product Data Sheet

Anti-Canine Distemper Virus Antibody (Capture)

Catalog #	Source	Reactivity	Applications
CAP1048-C	Mouse	CDV	ELISA, ICA, CLIA
Description	Mouse monoclonal capture antibody to Canine Distemper Virus		
Immunogen	Recombinant protein of Canine Distemper Virus		
Purification	The antibody was purified by immunogen affinity chromatography.		
Specificity	Recognizes native and recombinant Canine Distemper Virus protein.		
Clonality	Monoclonal		
Conjugation			
Form	IgG1, liquid in PBS, pH 7.4, 0.05% Proclin 300.		
Dilution			
Gene Symbol			
Alternative Names	CDV; Canine Distemper Virus		
Entrez Gene			
SwissProt			
Directions for Use	Recommended pair for sandwich immunoassay - Detection: CAP1048-D		
Storage/Stability	Shipped at 4°C. Upon delivery aliquot and store at -20°C for 2 years. Avoid freeze/thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com