

Product Data Sheet

RPL10 Blocking Peptide

| Catalog # | Source | Reactivity | Applications |
|---------------------------|---|-------------------------|--------------|
| CBP2016 | Synthetic | H, M, R, B, C, Mk, P, S | BL |
| Description | The peptide is used to block Anti-RPL10 Antibody (#CPA2016) reactivity. | | |
| Form | Lyophilized powder | | |
| Gene Symbol | RPL10 | | |
| Alternative Names | DXS648E; QM; 60S ribosomal protein L10; Laminin receptor homolog; Protein QM; Tumor suppressor QM | | |
| Entrez Gene | 6134 (Human); 110954 (Mouse); 81764 (Rat) | | |
| SwissProt | P27635 (Human); Q6ZWV3 (Mouse); Q6PDV7 (Rat) | | |
| Purity | >85% | | |
| Quality Control | The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. | | |
| Directions for Use | Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours. | | |
| Storage/Stability | Shipped at 4°C. Store at -20°C for one year. | | |

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com