

# Product Data Sheet

## HSP90 beta (pS254) Blocking Peptide

| Catalog #                 | Source  | Reactivity  | Applications |
|---------------------------|---|-------------|--------------|
| CBP3244                   | Synthetic   | H, M, R, Mk | BL           |
| <b>Description</b>        | The peptide is used to block Anti-HSP90 beta (pS254) Antibody (#CPA3244) reactivity.  |             |              |
| <b>Form</b>               | Lyophilized powder  |             |              |
| <b>Gene Symbol</b>        | HSP90AB1  |             |              |
| <b>Alternative Names</b>  | HSP90B; HSPC2; HSPCB; Heat shock protein HSP 90-beta; HSP 90; Heat shock 84 kDa; HSP 84; HSP84  |             |              |
| <b>Entrez Gene</b>        | 3326 (Human); 15516 (Mouse); 301252 (Rat)   |             |              |
| <b>SwissProt</b>          | P08238 (Human); P11499 (Mouse); P34058 (Rat)  |             |              |
| <b>Purity</b>             | >85%  |             |              |
| <b>Quality Control</b>    | The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.   |             |              |
| <b>Directions for Use</b> | Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours. |             |              |
| <b>Storage/Stability</b>  | Shipped at 4°C. Store at -20°C for one year.  |             |              |

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

### COHESION BIOSCIENCES LIMITED

**WEB**  
www.cohesionbio.com

**ORDER**  
order@cohesionbio.com

**SUPPORT**  
techsupport@cohesionbio.com

**CUSTOM**  
custom@cohesionbio.com