

Product Data Sheet

TECR Blocking Peptide

Catalog #	Source	Reactivity	Applications
CBP7228	Synthetic	H, M, R, B	BL
Description	The peptide is used to block Anti-TECR Antibody (#CPA7228) reactivity.		
Form	Lyophilized powder		
Gene Symbol	TECR		
Alternative Names	GPSN2; SC2; Very-long-chain enoyl-CoA reductase; Synaptic glycoprotein SC2; Trans-2,3-enoyl-CoA reductase; TER		
Entrez Gene	9524 (Human); 106529 (Mouse); 191576 (Rat)		
SwissProt	Q9NZ01 (Human); Q9CY27 (Mouse); Q64232 (Rat)		
Purity	>85%		
Quality Control	The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.		
Directions for Use	Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.		
Storage/Stability	Shipped at 4°C. Store at -20°C for one year.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com