


Angiotensinogen (1-14) (Porcine)

Catalog #	Source	Reactivity	Applications
CCP1946	Synthetic		
Description	Peptide to Angiotensinogen (1-14) (Porcine)		
Form	Lyophilized powder		
CAS Number	64315-16-8		
Molecular Formula	C ₈₅ H ₁₂₃ N ₂₁ O ₂₀		
Molecular Weight	1759.10		
Purity	> 95%		
Chemical Structure	Asp - Arg - Val - Tyr - Ile - His - Pro - Phe - His - Leu - Leu - Val - Tyr - Ser		
Storage/Stability	Shipped at 4°C. Store at -20°C for one year.		

Structure



Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com