

Product Data Sheet

ARF6 (2-13)

| Catalog # | Source | Reactivity | Applications |
|---------------------------|-----------|-----------------------------------------------------------------------|--------------|
| CCP1972 | Synthetic | | |
| Description | | Peptide to ARF6 (2-13) | |
| Form | | Lyophilized powder | |
| CAS Number | | | |
| Molecular Formula | | C60H102N16O17 | |
| Molecular Weight | | 1319.58 | |
| Purity | | > 95% | |
| Chemical Structure | | Gly - Lys - Val - Leu - Ser - Lys - Ile - Phe - Gly - Asn - Lys - Glu | |
| Storage/Stability | | Shipped at 4°C. Store at -20°C for one year. | |
| Structure | | | |

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com