

## Product Data Sheet

## Deslorelin Acetate

Catalog #	Source	Reactivity	Applications
CCP2249	Synthetic		
Description	Peptide to Deslorelin Acetate		
Form	Lyophilized powder		
CAS Number	57773-65-6		
Molecular Formula	C <sub>64</sub> H <sub>83</sub> N <sub>17</sub> O <sub>12</sub>		
Molecular Weight	1282.47		
Purity	> 95%		
Chemical Structure	Pyr - His - Trp - Ser - Tyr - D - Trp - Leu - Arg - Pro - NH <sub>2</sub>		
Storage/Stability	Shipped at 4°C. Store at -20°C for one year.		
Structure			

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

## COHESION BIOSCIENCES LIMITED

WEB  
[www.cohesionbio.com](http://www.cohesionbio.com)

ORDER  
[order@cohesionbio.com](mailto:order@cohesionbio.com)

SUPPORT  
[techsupport@cohesionbio.com](mailto:techsupport@cohesionbio.com)

CUSTOM  
[custom@cohesionbio.com](mailto:custom@cohesionbio.com)