

## Product Data Sheet

## Metamorphosin A

Catalog #	Source	Reactivity	Applications
CCP2619	Synthetic		
Description	Peptide to Metamorphosin A		
Form	Lyophilized powder		
CAS Number			
Molecular Formula	C34H47N9O8		
Molecular Weight	709.80		
Purity	> 95%		
Chemical Structure	Pyr - Gln - Pro - Gly - Leu - Trp - NH <sub>2</sub>		
Storage/Stability	Shipped at 4°C. Store at -20°C for one year.		
Structure			

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

## COHESION BIOSCIENCES LIMITED

WEB  
[www.cohesionbio.com](http://www.cohesionbio.com)

ORDER  
[order@cohesionbio.com](mailto:order@cohesionbio.com)

SUPPORT  
[techsupport@cohesionbio.com](mailto:techsupport@cohesionbio.com)

CUSTOM  
[custom@cohesionbio.com](mailto:custom@cohesionbio.com)