

Product Data Sheet

Preprogalanin 28-67 (Rat)

Catalog #	Source	Reactivity	Applications
CCP2803	Synthetic		
Description	Peptide to Preprogalanin 28-67 (Rat)		
Form	Lyophilized powder		
CAS Number			
Molecular Formula	C198H312N62O58		
Molecular Weight	4489.08		
Purity	> 95%		
Chemical Structure	Thr - Lys - Glu - Lys - Arg - Gly - Trp - Thr - Leu - Asn - Ser - Ala - Gly - Tyr - Leu - Leu - Gly - Pro - His - Ala - Ile - Asp - Asn - His - Arg - Ser - Phe - Ser - Asp - Lys - His - Gly - Leu - Thr - Gly - Lys - Arg - Glu - Leu - Pro		
Storage/Stability	Shipped at 4°C. Store at -20°C for one year.		
Structure			

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com