


Ser-Ala-SAP-IIB

Catalog #	Source	Reactivity	Applications
CCP2882	Synthetic		

Description	Peptide to Ser-Ala-SAP-IIB
Form	Lyophilized powder
CAS Number	140653-27-6
Molecular Formula	C ₄₂ H ₇₁ N ₁₃ O ₁₄ S ₂
Molecular Weight	1046.24
Purity	> 95%
Chemical Structure	Ser - Ala - Lys - Leu - Cys - Pro - Gly - Gly - Asn - Cys - Val (Disulfide bridge Cys5 - Cys10)

Storage/Stability Shipped at 4°C. Store at -20°C for one year.

Structure



Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com