

Biotin-Intermedin (Rat)

| Catalog # | Source | Reactivity | Applications |
|---------------------------|--|------------|--------------|
| CCP3020 | Synthetic | | |
| Description | Peptide to Biotin-Intermedin (Rat) | | |
| Form | Lyophilized powder | | |
| CAS Number | | | |
| Molecular Formula | C236H377N77O66S3 | | |
| Molecular Weight | 5445.29 | | |
| Purity | > 95% | | |
| Chemical Structure | Biotin - Pro - His - Ala - Gln - Leu - Leu - Arg - Val - Gly - Cys - Val - Leu - Gly - Thr - Cys - Gln - Val - Gln - Asn - Leu - Ser - His - Arg - Leu - Trp - Gln - Leu - Val - Arg - Pro - Ser - Gly - Arg - Arg - Asp - Ser - Ala - Pro - Val - Asp - Pro - | | |
| Storage/Stability | Shipped at 4°C. Store at -20°C for one year. | | |
| Structure | | | |

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, CHIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com