

Product Data Sheet

Parathyroid Hormone (39-84) (Human)

Catalog #	Source	Reactivity	Applications
CCP3073	Synthetic		
Description	Peptide to Parathyroid Hormone (39-84) (Human)		
Form	Lyophilized powder		
CAS Number	90880-43-6		
Molecular Formula	C211H357N67O72		
Molecular Weight	4984.5		
Purity	> 95%		
Chemical Structure	Ala - Pro - Leu - Ala - Pro - Arg - Asp - Ala - Gly - Ser - Gln - Arg - Pro - Arg - Lys - Lys - Glu - Asp - Asn - Val - Leu - Val - Glu - Ser - His - Glu - Lys - Ser - Leu - Gly - Glu - Ala - Asp - Lys - Ala - Asp - Val - Asn - Val - Leu - Thr - Lys - Ala - Lys - Ser - Gln		
Storage/Stability	Shipped at 4°C. Store at -20°C for one year.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com