

miRNA Mimic Negative Control-FAM

Catalog #	Source	Reactivity	Applications
CMH0000-FAM	Synthetic	N/A	
Description	Synthetic miRNA negative control-FAM for experiments using targeted miRNA transfection		
Specificity	miRNA Mimics are chemically modified, double-stranded miRNA-like RNA which are designed to copy the functionality of mature endogenous miRNA upon transfection. At the 5'-end, it is synthesized with a partially complementary motif to 3'UTR end of the target gene, which allows the miRNA mimic to specifically bind to the target. Transfection of mimics followed by downstream gene expression analysis or phenotypic analysis, is performed to elucidate the targets and roles of particular miRNAs.		
Form	Lyophilized powder		
Gene Symbol			
Accession No.			
Components	This synthetic miRNA is based on the mature miRNA sequence conjugating with FAM. It does not contain the full precursor miRNA stem-loop.		
Directions for Use	We recommend re-suspending the lyophilized synthetic miRNA using DNase and RNase-free ddH ₂ O. To make a 100 uM stock, dissolve the lyophilized powder using 50 ul of ddH ₂ O.		
Storage/Stability	Shipped at 4 °C. Store at -20 °C for one year. Avoid freeze-thaw cycles after reconstitution.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED**WEB**www.cohesionbio.com**ORDER**order@cohesionbio.com**SUPPORT**techsupport@cohesionbio.com**CUSTOM**custom@cohesionbio.com