

Product Data Sheet

Anti-Hepatitis B Virus Core Antigen Antibody-FITC labeled

Catalog #	Source	Reactivity	Applications
CPA9940	Mouse	Hepatitis B virus	E
Description	Mouse monoclonal antibody to Hepatitis B Virus Core Antigen-FITC labeled		
Immunogen	Recombinant Hepatitis B Virus Core Antigen		
Purification	Immunogen affinity purified		
Specificity	Reacts only with Hepatitis B Virus Core Antigen.		
Clonality	Monoclonal		
Conjugation	FITC		
Form	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.05% Sodium Azide.		
Dilution	ELISA: Use at an assay dependent dilution.		
Gene Symbol			
Alternative Names	Capsid protein; Core and e antigen; Core antigen; Core antigen; Core protein; HBc; HBcAg; HBVgp4; Hepatitis B Virus core antigen; p21.5; precore/core protein		
Entrez Gene			
SwissProt			
Storage/Stability	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com