

# Product Data Sheet

## Anti-Casein Antibody-Biotin labeled

Catalog #	Source	Reactivity	Applications
CPA9956	Rabbit	N/A	E
<b>Description</b>	Rabbit polyclonal antibody to Casein-Biotin labeled		
<b>Immunogen</b>	Full length Casein		
<b>Purification</b>	Immunogen affinity purified		
<b>Specificity</b>	Reacts only with Casein.		
<b>Clonality</b>	Polyclonal		
<b>Conjugation</b>	Biotin		
<b>Form</b>	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.05% Sodium Azide.		
<b>Dilution</b>	ELISA: Use at an assay dependent dilution.		
<b>Gene Symbol</b>			
<b>Alternative Names</b>			
<b>Entrez Gene</b>			
<b>SwissProt</b>			
<b>Storage/Stability</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.		

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

### COHESION BIOSCIENCES LIMITED

**WEB**  
www.cohesionbio.com

**ORDER**  
order@cohesionbio.com

**SUPPORT**  
techsupport@cohesionbio.com

**CUSTOM**  
custom@cohesionbio.com