

Product Data Sheet

Anti-Human Serum Albumin Antibody-Biotin labeled

Catalog #	Source	Reactivity	Applications
CPA9958	Rabbit	H	E
Description	Rabbit polyclonal antibody to Human Serum Albumin-Biotin labeled		
Immunogen	Human Serum Albumin		
Purification	Immunogen affinity purified		
Specificity	Reacts only with Human Serum Albumin.		
Clonality	Polyclonal		
Conjugation	Biotin		
Form	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, 50% glycerol, 0.05% Sodium Azide.		
Dilution	ELISA: Use at an assay dependent dilution.		
Gene Symbol	ALB		
Alternative Names	Serum albumin		
Entrez Gene	213 (Human)		
SwissProt	P02768 (Human)		
Storage/Stability	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com