

# Product Data Sheet

## TMB ELISA Substrate

Catalog #	Source	Reactivity	Applications
CRG1113			ELISA
Description	TMB ELISA Substrate detects horseradish peroxidase (HRP) activity.		
Form	Liquid		
Directions for Use	<p>(a) After incubation with horseradish peroxidase (HRP) conjugated secondary antibody, wash plates with PBS/0.05% Tween-20.</p> <p>(b) Add 100 µl of TMB solution/well. Incubate for a specific amount of time or until desired signal is achieved.</p> <p>(c) Add 100 µl of ELISA Stop Solution.</p> <p>(d) Read at Absorbance 450 nm.</p>		
Components	TMB is a one component system that is ready to use from the bottle.		
Storage/Stability	Store at 4 °C in dark for 1 year, do not freeze.		

### COHESION BIOSCIENCES LIMITED