

# Product Data Sheet

## Recombinant Mouse IGHG1 Protein

Catalog #	Source	Reactivity	Applications
CRP1726	Human cells	Mouse	E, WB, SDS-PAGE, MS
<b>Description</b>	Recombinant Mouse IGHG1 Protein is produced by our mammalian expression system and the target gene encoding Pro99-Lys324 is expressed.		
<b>Form</b>	Lyophilized from a 0.2 µM filtered solution of 20mM PB, 150mM NaCl, pH 7.4.		
<b>Gene Symbol</b>	IGHG1		
<b>Alternative Names</b>	Ig gamma-1 chain C region		
<b>Entrez Gene</b>	16017 (Mouse)		
<b>SwissProt</b>	P01868 (Mouse)		
<b>Purity</b>	Greater than 95% as determined by reducing SDS-PAGE.		
<b>Chemical Structure</b>	GSGSPRDCGC KPCICTVPEV SSVFIFPPKP KDVLITLTP KVTCVVVDIS KDDPEVQFSW FVDDDEVHTA QTQPREEQFN STFRSVSELP IMHQDWLNGK EFKCRVNSAA FPAPIEKTIS KTKGRPKAPQ VYTIPPPKEQ MAKDKVSLTC MITDFFPEDI TVEWQWNGQP AENYKNTQPI MDTDGSYFVY SKLNVQKSNW EAGNTFTCSV LHEGLHNHHT EKSLSHSPGK		
<b>Quality Control</b>	Endotoxin: Less than 0.1 ng/µg (1 IEU/µg) as determined by LAL test.		
<b>Directions for Use</b>	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 µg/ml. Dissolve the lyophilized protein in 1X PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
<b>Storage/Stability</b>	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 2-8°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.		

**Application key:** E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

**Species reactivity key:** H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

### COHESION BIOSCIENCES LIMITED

**WEB**  
www.cohesionbio.com

**ORDER**  
order@cohesionbio.com

**SUPPORT**  
techsupport@cohesionbio.com

**CUSTOM**  
custom@cohesionbio.com