

Product Data Sheet

Glycerol-3-Phosphate Oxidase

Catalog #	Source	Reactivity	Applications
CZM1034	microorganism		
Description	Recombinant Glycerol-3-Phosphate Oxidase from microorganism		
Specificity	>15 U/mg		
Form	Lyophilized powder		
Alternative Names	EC 1.1.3.21; GPO; glycerol phosphate oxidase; glycerol-1-phosphate oxidase; glycerol phosphate oxidase; L-a-glycerophosphate oxidase; a-glycerophosphate oxidase; L-a-glycerol-3-phosphate oxidase		
CAS Number	9046-28-0		
Molecular Weight			
Purity	>90%		
Unit Definition			
Storage/Stability	Shipped at 4° C. Store at -20° C for 2 years.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com