

Product Data Sheet

Cystathionine beta-Synthase

Catalog #	Source	Reactivity	Applications
CZM1082	E. coli		
Description	Cystathionine beta-Synthase from E. coli		
Specificity	>40 U/mg protein, >5 U/mg powder		
Form	Light yellow lyophilized powder		
Alternative Names	EC 4.2.1.22; CBS; beta-thionase; methylocysteine synthase; serine sulfhydrase		
CAS Number	9023-99-8		
Molecular Weight			
Purity	>90%		
Storage/Stability	Shipped at 4°C. Store at -20°C for 2 years.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com